



## Recommendations

RECOMMENDATIONS	ACTIONS
For the 1 <sup>st</sup> few weeks of using	Hover on Dragon bar
Dragon, to ensure Dragon learns	Click Audio
and remembers your verbal	Click Launch Accuracy Tuning
patterns, activate the Launch	Acoustic and Language Model Optimizer screen
Accuracy Tuning <u>weekly</u>	displayed
	Click <b>Go</b>
After a few weeks, activate the	Follow screen prompts
Launch Accuracy Tuning on a	
<u>monthly</u> basis	
Be and this illustration of a state	
<b>Beware</b> – this will slow down the	
whole laptop for approximately 30-60mins	
Do a Monthly Backup	Hover on Dragon bar
· · · · / · · · ·	Click <b>Profile</b> icon
	Click Backup User Profile
	Wait
	When the Backup is complete, the Dragon bar will
	return to normal i.e. no green message displayed
To Improve Dictation Accuracy	You can read various dictation samples, to help
	Dragon recognise & improve the quality of its memory
	of your voice
	Hover on Dragon bar
	Click <b>Audio</b> icon
	Click Read text to improve accuracy
	Training Wizard window displayed
	Click Next
	Click on required Dictation sample
	Click Next